



Operation mat for footstool and operating room floor – Product Data Sheet

Provider:	Svenska Good Medical AB
Phone number:	+46 31 785 19 61
Product name:	AireRx Op Mat S
Article number:	AR001
Area of usage:	Pressure relieving operation mat for footstool
Material:	EV50: water absorption <1 %, recommended maximum operating temperature 65 °C. For more information please visit: http://www.zotefoams.com/wp-content/uploads/2015/12/EV50-Dec-2015.pdf
Measures:	445 x 345 mm
Disposal:	After use, handled consumed article in accordance with the hospital instructions.
	Do Not Reuse / Single Use Only

Operationsmatta för pall och operationssalsgolv – Produktdatablad

Leverantör:	Svenska Good Medical AB
Telefonnummer:	+46 31 785 19 61
Produktnamn:	AireRx Op Mat S
Artikelnummer:	AR001
Användningsområde:	Tryckavlastande matta för operationspall
Material:	EV50: ej vattenabsorberande (<1 %), tål temperatur upp till 65 °C. För mer information var god se: http://www.zotefoams.com/wp-content/uploads/2015/12/EV50-Dec-2015.pdf
Mått:	445 x 345 mm
Kassering:	Efter användning hanteras förbrukad artikel i enlighet med sjukhusets destrueringsanvisningar.
	Endast för engångsbruk.